



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Auto Urine Analyzer MA-240

FEATURES:

- High precision sampler
- Auto mixing before aspirating samples
- Particular IC card control system and system upgrade function
- STAT: emergency insert dedicated location

TESTING ITEMS:

Urobilinogen (UBG), Bilirubin (BIL), Ketone body (KET), Occult blood (BLD), Protein (PRO), Nitrite (NIT), Leukocyte (LEU), Glucose (GLU), Specific Gravity (SG), Power of Hydrogen (pH), Vitamin C (VC), Micro albumin (MAL), Creatinine (CRF), Calcium ion (CAL)



MODEL SPECIFICATIONS:

Model	MA-240
Assay Method	Albedometry
Test Speed	240 tests/hr
Sample Rack	5 test tube racks (5*10 samples)
Sample Volume	Not lower than the marked scale line on the test tube rack
Sampling Volume	<1.0ml
Wavelength	525nm, 610nm, 660nm
Storage	20 000 results
Interface	RS-232
Language	English and Spanish. Other languages on request
Printer	Built-in thermal printer
Working Environment	Temperature: 15~35°C, humidity: ≤75%
Power Supply	AC220V±10% 50/60Hz; 110V±10% 60Hz
Fuse Specification	250V, 2A
Package Size (W*D*H)	800*595*995mm
Gross Weight (kg)	83