



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Automatic Uncovering Centrifuge TLS-AU4 / TLS-AU5 / TLS-AU5R

FEATURES:

- LCD display, microprocessor control
- Stainless steel chamber
- Mechatronics motor door lock
- Can store 20 procedures
- 10 levels speed increase and decrease control
- Brushless motor: variable frequency, free from contamination caused by carbon dust



TLS-AU4



TLS-AU5/TLS-AU5R



48*7/5ml



76*7/5ml



124*7/5ml

MODEL SPECIFICATIONS:

Model	TLS-AU4	TLS-AU5	TLS-AU5R
Capacity	48*2/5ml; 76*2/5ml	124*2/5ml; 148*5ml	124*2/5ml; 148*5ml
Max. Speed	4000rpm	4000rpm	4000rpm
Max. RCF	3100×g	3400×g	3400×g
Speed Precision	±50rpm		
Refrigeration System	/		CFC-free (R404a)
Temp. Range	/		-20°C~40°C
Temp. Precision	/		±1°C
Control and Drive System	Frequency motor, micro control, direct driver		
Identify Rotor	/		Automatic
Timing Range	0~99h59min		
Noise	≤65dB		
Consumption	800W	1000W	2000W
Power Supply	AC110/220V, 50/60Hz		
External Size (W*D*H)	630*500*420mm	710*560*770mm	710*560*770mm
Package Size (W*D*H)	780*840*570mm	900*750*970mm	900*750*1195mm
Nett / Gross Weight (kg)	50 / 88	75 / 105	150 / 215