



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Blood Centrifuge TLS-TB12

#### FEATURES:

- Suitable for blood group serology, blood type routine test, microsurgery Gel, erythrocyte washing and immunoassay, etc



Lymphocyte Cleaning  
HLS Rotor



Blood Type  
Card Rotor



Erythrocyte Washing  
SERO Rotor



**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>TLS-TB12</b>
Max. Speed	4000rpm
Max. RCF	2200×g
Consumption	200W
Power Supply	110/220V, 50/60Hz
External Size (W*D*H)	500*450*280mm
Package Size (W*D*H)	700*600*450mm
Net / Gross Weight (kg)	35 / 40

**Rotor:**

<b>Rotor</b>	<b>Centrifugal Time</b>	<b>RCF (xg)</b>	<b>Applications</b>
Erythrocyte Washing SERO Rotor	60min	500	Blood type determination, hemagglutination reflect observed
	15min	1000	Cross-suitability experiment, engage ball prime experiment
	60min	1000	Blood cell washing, extraction of serum and plasma
Lymphocyte Cleaning HLS Rotor	180min	2000	Lymphocyte separation, cultured cell separation
	3min	1000	Removal of platelets (thrombin treatment)
	60min	1000	Lymphocyte cleaning
Blood Type Card Rotor	/	2500	6 and 8 associated micro blood type card, 24 cards, centrifuge holds 24 ScanGel cards