

CENTRIFUGES

CRUDE OIL CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Crude Oil Centrifuge MA-OIL5B

#### FEATURES:

- 3-phase centrifuge
- Stainless steel chamber
- Electrical lid locking (over-speed and over-temperature protection)
- Programmable microprocessor control system
- LCD display shows rotor number, speed, time, temperature, RCF, rate of speed as well as error information
- Adaptable band heater (rapid heating speed, uniform temperature, temperature-control and constant function)
- Brushless motor



4\*100ml

MA-OIL5B-03

**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>MA-OIL5B</b>
Max. Speed	4000rpm
Max. RCF	2498xg
Max. Capacity	4*100ml
Timer Range	1~99min
Noise	≤60dB
Power Supply	AC110/220V, 50/60Hz
External Size (W*D*H)	685*500*405mm
Package Size (W*D*H)	870*590*520mm
Net Weight (kg)	67
Gross Weight (kg)	80