



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

DispensMate Bottle-Top Dispenser DM-004B

FEATURES:

- Chemical resistance
- Fully Autoclavable at 121°C
- Cover a volume range from 0.5ml to 50ml
- Easy to clean and maintain
- Optional flexible discharge tube with safety handle
- Made of PTFE, FEP, BSG, PP
- Vapour pressure max. 500mbar, viscosity max. 500mm²/s, temperature max. 40°C, density max. 2.2g/cm³
- DispensMate plus is supplied with S40, GL32, GL38, GL25, GL28 sized adapters



MODEL SPECIFICATIONS:

Model	DM-004B				
Volume Range	Graduation	A \leq		CV \leq	
		%	μ l	%	μ l
0.5-5ml	0.1ml	0.5	25	0.1	5
1.0-10ml	0.2ml	0.5	50	0.1	10
2.5-25ml	0.5ml	0.5	125	0.1	25
5.0-50ml	1.0ml	0.5	250	0.1	50

A = Accuracy; CV = Coefficient of variation

ORDERING INFORMATION:

Model No	Range Volume
DM004B5	0.5-5ml
DM004B10	1-10ml
DM004B25	2.5-25ml
DM004B50	5-50ml

Model No	Volume Range
Standard Accessories	
DM004BA1	S40, Adapter, 45/40mm
DM004BA2	GL32, Adapter, 45/32mm
DM004BA3	GL38, Adapter, 45/38mm
DM004BA4	GL25, Adapter, 32/25mm
DM004BA5	GL28, Adapter, 32/28mm
DM004BFT	Filling tube (m)
DM004BDT	Discharge tube (m)
Optional Accessories	
DM004BRB	Reagent Bottle (Brown, 1L)