



+27 60 999 9140
 info@microanalytica.co.za
 489 Jacqueline Drive, Garsfontein
 Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Electronic Bottle-Top Dispenser

FEATURES:

- Motor operation minimizes the grip force
- Electronic control decreases the risk of repetitive strain injury
- Remote control panel prevents manual intervention during operation
- Two dispensing modes
- Dispenser
- Stepper function



MODEL SPECIFICATIONS:

Model	EB004
Volume Range	0.01~99.99ml, max. piston lift is 10ml, resolution is 10µl
Volume Accuracy	R =0.2%; CV = 0.07%
Operating Temperature Range	10°C~30°C
Quality Control	DIN EN ISO 8655
Control Type	External control, including the stirrer and the burette program