



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Floor Standing Large Capacity High Speed Refrigerated Centrifuge FHS-VH10RL

FEATURES:

- LCD display, microprocessor control, indicating setting parameters and actual parameters at the same time
- AV motor: high torque and variable frequency, free from contamination caused by carbon dust
- Imported drive system and automatic control system to ensure steady and reliable operation
- Safety: Over-temperature and over-speed protection, electronic locking, imbalance detection
- Can set 16 procedures, parameters can be modified during running
- 10 levels increase and decrease speed control
- Can indicate rotor speed ascending and declining curves, centrifugal integral curves and temperature curves



6*1000ml

MODEL SPECIFICATIONS:

Model	FHS-VH10RL
Max. Capacity	6*1000ml
Max. Speed	10000rpm
Max. RCF	18300xg
Speed Precision	±50rpm
Control and Driver System	Frequency motor, micro control, direct driver
Refrigeration System	CFC-free (R404a)
Temp. Range	-20°C~40°C
Temp. Precision	±1°C
Timing Range	0~99h59min
Noise	<65dB
Consumption	5.7KW
Power Supply	AC110/220V, 50/60Hz
External Size (W*D*H)	840*730*950mm
Package Size (W*D*H)	1085*835*1165mm
Net / Gross Weight (kg)	280 / 370

ROTOR:

Rotor	Capacity	Speed (rpm)	RCF (xg)
Angle Rotor	6*1000ml	7000	11500
	6*300ml	10000	18300
Swing Rotor	6*500ml	8000	11800
	6*1000ml	4200	5100