

CENTRIFUGES

GEL CARD CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Gel Card Centrifuge MA-GC12 / MA-GC24

FEATURES:

- Gel card centrifuge is used in blood serology, blood type routine test, erythrocyte washing, micro column agglutination, immunology and other tests
- DC Brushless motor, low noise
- Microprocessor control system & soft touch control panel
- LCD display: RCF, Speed, Time
- Electrical lid interlocking, over-speed and imbalance protection
- 20 standard and programmable programs are available



12 cards



24 cards



MODEL SPECIFICATIONS:

Model	MA-GC12	MA-GC24
Max. Speed	4000rpm	
Max. RCF	2250xg	
Rotor	12 cards	24 cards
Speed Precision	±20rpm	
Timing Range	0~99min	
Noise	≤55dB	
Power Supply	AC110/220V±10%, 50/60Hz	
External Size (W*D*H)	460*320*230mm	500*360*250mm
Package Size (W*D*H)	560*395*330mm	610*435*340mm
Gross Weight (kg)	24	30