

HEATING BLOCKS



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Heating Block MA-III

FEATURES:

- Wide range of interchangeable aluminium blocks
- Heat to 105°C, optional external temperature for all applications - tubes, vials and microplates
- Blocks optional available as accessories for all applications - tubes, vials and microplates
- Block extraction tool is supplied, allowing blocks to be removed easily



MODEL SPECIFICATIONS:

Model	MA-III
Temperature Control Range	RT+5°C ~ 100°C
Time required for temperature to increase from 20°C to 100°C	≤15min
Temperature Control Discrepancy	≤±0.5°C
Display Discrepancy	±0.1°C
Homogeneity of Block Temperature	≤± 0.5°C
Maximum Heating Power	100W
Maximum Temp.	100°C
Time Range	1min~99h59min
External Size (mm)	225 *160 *197
Package Size (mm)	330*280*340
Net Weight (kg)	3
Gross Weight (kg)	6