

CENTRIFUGES

HAEMATOCRIT CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Haematocrit Centrifuge HEM-HC12B

FEATURES:

- Micro brushless motor, maintenance, no powder pollution, micro-computer control, LCD displays time and speed
- Electrical lid locking
- Over-speed and imbalance protection
- Frame is 3 tiers steel jacket
- Stainless steel chamber



MA-HEM/HC12B-03

MODEL SPECIFICATIONS:

Model	HEM-HC12B
Max. Speed	12000rpm
Max. RCF	15800xg
Max. Capacity	24 capillaries
Timer Range	0~99min59s
Noise	≤55dB
Power Supply	AC110/220V, 50/60Hz
External Size (W*D*H)	310*270*220mm
Package Size (W*D*H)	395*310*260mm
Nett Weight (kg)	14
Gross Weight (kg)	16

ROTOR:

Rotor No.	Rotor	Capacity	Speed (rpm)	RCF (xg)
No. 31001	Haematocrit rotor	24 capillaries	12000	15800