

## PIPETTES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Levo Plus Pipette Filler

#### FEATURES:

- Volume range 0.1-100ml
- One hand operation
- Lightweight, ergonomic design
- Efficient lithium-ion battery
- Low battery alarm
- Powerful motor fills a 25ml pipette within 5 seconds
- Integral and replaceable 0.45µm filter
- Large LCD display
- Eight speeds are available to aspirate and dispense liquid



**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>LP7033110100</b>
Aspiration Speeds	25ml<5s (6 shift)
Dispensing Speed	Motor (6 shift) / Gravity
Battery	Lithium-ion
Battery Run Time	8 hours of intermittent use
Charging Time	Approx. 2 hours
Pipette	Glass or Plastic Pipette (0.1-100ml), Pasteur Pipettes
Filter	0.45µm Hydrophobic

**ORDERING INFORMATION:**

<b>Model No</b>	<b>Description</b>
LP7033110100	Levo Plus with AC adapter, spare 0.45µm filter and wall stand, Euro plug
	Levo Plus with AC adapter, spare 0.45µm filter and wall stand, North America plug
	Levo Plus with AC adapter, spare 0.45µm filter and wall stand, UK plug
<b>Standard Accessories</b>	
LP7033110100C	Levo Plus charger
LP7033110100B	Levo Plus Lithium battery
LP7033110100F	Filter 0.45µm, 1pc/pk