



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Micro Table Top High Speed Centrifuge THS-MH16

FEATURES:

- LCD display: indicating setting parameters and actual parameters, including RCF
- Variable frequency brushless motor
- Stainless steel chamber
- Mute mechatronics motor door lock
- Can store 20 procedures
- 10 levels speed increase and decrease control



12*1.5/2.0ml



24*1.5/2.0ml



10*5ml



12*10ml



MODEL SPECIFICATIONS:

Model	THS-MH16
Max. Capacity	12*10ml
Max. Speed	16500rpm
Max. RCF	18780×g
Timing Range	0~99h59min
Control and Drive System	Frequency motor, micro control, direct drive
Noise	≤60dB
Consumption	350W
Power Supply	AC110/220V, 50/60Hz
External Size (W*D*H)	380*260*310mm
Package Size (W*D*H)	620*360*410mm
Net / Gross Weight (kg)	20 / 28

ROTOR:

Rotor	Capacity	Speed (rpm)	RCF (×g)
Angle Rotor	12*1.5/2.0ml	16500	18780
	18*1.5/2.0ml	15000	17780
	24*1.5/2.0ml	14000	18180
	10*5ml	14000	13520
	8*7ml	14000	12340
	12*10ml	12000	13400