

PRECISION BALANCES

PB-DRP ROUND PAN MODEL RANGE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

PB-DRP Series Round Pan Electronic Balance (0.1g/0.01g)

FEATURES:

- High precision load cell
- Aluminum alloy die cast base and mains
- Stainless steel pan
- LCD display
- Full capacity subtraction
- Stabilization time 1.5 seconds
- Overload protection
- Below balance weighing facility
- Checking weighing, piece counting function
- RS232 Interface



MODEL SPECIFICATIONS:

Model No	PB2001DRP	PB6001DRP	PB8001DRP	PB10001DRP	PB12001DRP	PB10002DRP	PB20002DRP
Capacity	0-200g	0-600g	0-800g	0-1000g	0-1200g	0-1000g	0-2000g
Readability	0.1g/0.01g						
Repeatability	±0.2g/±0.02g						
Linearity	±0.2g/±0.02g						
Output Interface	RS232C						
Pan Size	Φ125mm						
Dimensions (L*W*H)	29cmx18.5cmx8cm						
Packing Size (L*W*H)	30.2cmx22cmx14.5cm						
Net weight	1800g						
Gross weight	2200g						
Power supply	110V/220-240V AC						