

## PRECISION BALANCES

### PB-DSP SQUARE PAN MODEL RANGE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

#### PB-DSP Series Square Pan Electronic Balance (1g/0.1g/0.01g)

##### FEATURES:

- High precision load cell
- Aluminum alloy die cast base and mains
- Stainless steel pan
- LCD display
- Full capacity subtraction
- Stabilization time 1.5 seconds
- Overload protection
- Below balance weighing facility
- Checking weighing, piece counting function
- RS232 Interface



**MODEL SPECIFICATIONS:**

<b>Model No</b>	<b>PB2000DSP PB20001DSP</b>	<b>PB3000DSP PB30001DSP PB30002DSP</b>	<b>PB4000DSP PB40001DSP PB40002DSP</b>	<b>PB5000DSP PB50001DSP PB50002DSP</b>	<b>PB6000DSP PB60001DSP</b>	<b>PB10000DSP PB100001DSP</b>
Capacity	0-2000g	0-3000g	0-4000g	0-5000g	0-6000g	0-10000g
Readability	1g/0.1g/0.01g					
Repeatability	±1g/±0.2g/±0.02g					
Linearity	±1g/±0.2g/±0.02g					
Output Interface	RS232C					
Pan Size	170mmx170mm					
Dimensions (L*W*H)	29cmx18.5cmx8cm					
Packing Size (L*W*H)	30.2cmx22cmx14.5cm					
Net Weight	2200g					
Gross Weight	2600g					
Power Supply	110V/220-240V AC					