

## PRECISION BALANCES

### PB-DWS WINDSHIELD MODEL RANGE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

#### PB-DWS Series Windshield Electronic Balance (0.01g)

##### FEATURES:

- High precision load cell
- Aluminum alloy die cast base and mains
- Stainless steel pan
- LCD display
- Full capacity subtraction
- Stabilization time 1.5 seconds
- Overload protection
- Below balance weighing facility
- Checking weighing, piece counting function
- RS232 Interface
- Plastic windshield as standard



**MODEL SPECIFICATIONS:**

<b>Model No</b>	<b>PB1002DWS</b>	<b>PB2002DWS</b>	<b>PB3002DWS</b>	<b>PB5002DWS</b>	<b>PB6002DWS</b>
Capacity	0-100g	0-200g	0-300g	0-500g	0-600g
Readability	0.01g				
Repeatability	±0.02g				
Linearity	±0.02g				
Output Interface	RS232C				
Pan Size	Φ125mm				
Dimensions (L*W*H)	29cmx18.5cmx8cm				
Packing Size (L*W*H)	30.2cmx22cmx14.5cm				
Net weight	2000g				
Gross weight	2400g				
Power supply	110V/220-240V AC				