

PRECISION BALANCES

PB-Y MODEL RANGE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### PB-Y Series Electronic Balance (0.01g / 0.1g)

#### FEATURES:

- High precision load cell
- Rugged plastic housing
- Big size stainless steel pan
- Super bright LCD display with backlight
- Built-in rechargeable battery
- Height adjustable feet
- Full capacity subtraction
- Stabilization time 1.5 seconds
- Overload protection
- Below weighing facility
- Checking weighing, piece counting function
- Bubble level for choice



**MODEL SPECIFICATIONS:**

<b>Model No</b>	<b>PB1002Y</b>	<b>PB2002Y</b>	<b>PB3002Y</b>	<b>PB5002Y</b>	<b>PB6002Y</b>	<b>PB10002Y</b>	<b>PB20002Y</b>
Capacity	0-100g	0-200g	0-300g	0-500g	0-600g	0-1000g	0-2000g
Readability	0.01g						
Repeatability	±0.02g						
Linearity	±0.02g						
Pan Size	128mmx157mm						
<b>Model No</b>	<b>PB10001Y PB1000Y</b>	<b>PB15001Y PB1500Y</b>	<b>PB20001Y PB2000Y</b>	<b>PB30001Y PB3000Y</b>	<b>PB40001Y PB4000Y</b>	<b>PB50001Y PB5000Y</b>	<b>PB60001Y PB6000Y</b>
Capacity	0-1000g	0-1500g	0-2000g	0-3000g	0-4000g	0-5000g	0-6000g
Readability	0.1g/1g						
Repeatability	±0.2g/±1g						
Linearity	±0.2g/±1g						
Pan Size	143mmx192mm						
Output Interface	RS232C						
Power Supply	AC 110 – 240V Adapter						