

CENTRIFUGES

MINI-CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Single Speed Mini Centrifuge MA Mini-7 / MA Mini-5

Used for serum extraction from whole blood, extraction from different samples, micro blood cell separation, bio-sample treatment, samples quick sedimentation and PCR sample treatment.

#### FEATURES:

- Rotor with 8 positions and rotor for PCR strip tube
- Unique rotor clip design
- Low noise, stable running



**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>MA Mini-7</b>	<b>MA Mini-5</b>
Max. Speed	7000rpm	5000rpm
Max. RCF	2650xg	1350xg
Capacity	0.2/0.5/1.5/2.0ml x 8; 8 x 4 x 0.2ml (PCR tube)	
Rotor	8*2.0/1.5/0.5/0.2ml rotor (0.2ml and 05.ml adaptor), 8*4*0.2ml PCR tube	
Power Supply	AC110V – 265V, 50/60Hz	
Consumption	20W	
Noise	45dB(A)	
External Size (W*D*H)	178*178*115mm	
Package Size (W*D*H)	320*250*210mm	
Net Weight (kg)	1	
Gross Weight (kg)	2	