

CENTRIFUGES

TABLE TOP DAIRY CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table Top Dairy Centrifuge TTD-MF5A

FEATURES:

- Special equipment for milk and dairy analysis, such as analysis of the samples after fat separation through pasteurization, Gerber's method, etc
- The LCD display indicates the speed, time and temperature
- Microprocessor control
- Electrical lid locking, over-speed, over-temperature and imbalance protection
- Stainless steel centrifugal chamber and lid



TTD-MF5A-03

MODEL SPECIFICATIONS:

Model	TTD-MF5A
Max. Capacity	8*Gerber tubes
Max. Speed	1500rpm
Max. RCF	420xg
Speed Precision	±50rpm
Temp. Range	RT+10°C~65°C
Timing Range	0~99h59min
Noise	≤60dB
Consumption	1.5KW
Power Supply	AC110/220V, 50/60Hz
External Size (W*D*H)	520*415*330mm
Package Size (W*D*H)	720*610*530mm
Net / Gross Weight (kg)	35 / 60

Rotor:

No.	Rotor	Capacity	Speed (rpm)	RCF (xg)
No. 1	Angle Rotor	8*Gerber tubes	1500	420