

CENTRIFUGES

TABLE TOP HIGH SPEED CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table Top High Speed Centrifuge THS-TH16 / THS-TH18

FEATURES:

- LCD display: indicating settings parameters and actual parameters, including RCF
- Variable frequency brushless motor
- Stainless steel chamber
- Mute mechatronics motor door lock
- Can store 20 procedures



12*1.5/2.0ml



6*50ml



MODEL SPECIFICATIONS:

Model	THS-TH16	THS-TH18
Max. Capacity	4*100ml	
Max. Speed	16500rpm	18500rpm
Max. RCF	18780×g	23610×g
Timing Range	0~99h59min	
Noise	≤60dB	
Consumption	050W	
Power Supply	AC110/220V, 50/60Hz	
External Size (W*D*H)	500*440*330mm	710*630*350
Package Size (W*D*H)	750*640*530mm	
Net / Gross Weight (kg)	30 / 50	

ROTOR:

Rotor	THS-TH16			THS-TH18		
	Capacity	Speed (rpm)	RCF (xg)	Capacity	Speed (rpm)	RCF (xg)
Angle Rotor	12*1.5/2.0ml	16500	18780	12*1.5/2.0ml	18500	23610
	18*1.5/2.0ml	15000	17770	18*1.5/2.0ml	15000	17770
	24*1.5/2.0ml	14000	18180	24*1.5/2.0ml	14000	18180
	10*5ml	13500	13520	10*5ml	14000	13520
	8*15ml	12000	14800	8*15ml	12000	14800
	12*10ml	12000	13400	12*10ml	12000	13400
	6*50ml	11000	11610	6*50ml	11000	11610
	6*50ml	6000	3460	6*50ml	6000	3460
Microplate	2*2*48 Well	3000	1400	2*2*48 Well	3000	1400