



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table Top Low Speed Centrifuge TLS-TL4 / TLS-TL4II / TLS-TL4III

FEATURES:

- LCD display, microprocessor control
- Brushless motor: variable frequency
- Multiple shock absorber layers, automatic balance of rotor to achieve smooth running
- Can store 20 procedures



TLS-TL4



TLS-TL4II / TLS-TL4III



8*20ml



24*10ml

MODEL SPECIFICATIONS:

Model	TLS-TL4	TLS-TL4II	TLS-TL4III
Capacity	8*20ml	12*20ml, 18*10ml, 24*10ml	4*50ml, 12*15ml
Max. Speed	4000rpm		
Max. RCF	2200xg	2600xg	2100xg
Speed Precision	±30rpm		
Timing Range	0~99h59min inching function		
Noise	≤60dB		
Power Supply	AC110/220V, 50/60Hz		
External Size (W*D*H)	380*260*310mm	500*450*280mm	
Package Size (W*D*H)	490*610*510mm	700*630*480mm	
Net Weight (kg)	18	20	
Gross Weight (kg)	33	35	