

CENTRIFUGES

TABLE TOP LOW SPEED CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table Top Low Speed Centrifuge TLS-TL4M

FEATURES

- Stainless steel chamber
- Microprocessor control, brushless DC motor
- LCD display of running parameters
- Unbalance protection, automatic electrical lid locking



12*15ml



18*10ml



24*10ml

MODEL SPECIFICATIONS:

Model	TLS-TL4M
Max. Speed	4000r/min
Max. RCF	2200×g
Max. Capacity	24*10ml
Speed Precision	±20rpm
Timing Range	0~999min
Power Supply	AC 220V/110V, 50Hz/60Hz, 5A
Rated Power	250W
Noise	≤65dB
External Size (W*D*H)	350*420*270mm
Net Weight (kg)	18
Package Size (W*D*H)	410*460*330mm
Gross Weight (kg)	22

ROTOR:

No.	Rotor	Capacity	Max. Speed (r/min)	Max. RCF (xg)
No. 1	Angle Rotor	12*15ml	4000	2200
No. 2		18*10ml	4000	2200
No. 3		24*10 ml	4000	2200