

CENTRIFUGES

TABLE TOP LOW SPEED CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table Top Low Speed Centrifuge TLS-TL6M

FEATURES

- Stainless steel chamber
- Microprocessor control, brushless DC motor
- Display running parameters
- Unbalance protection, automatic electric lid locking



4*250ml



MA-TLS/TL6M-03

MODEL SPECIFICATIONS:

Model	TLS-TL6M
Max. Speed	6000rpm
Max. RCF	5310×g
Max. Capacity	4*500ml
Speed Precision	±30rpm
Time Range	0~999min
Power Supply	AC 220V/110V, 50/60Hz 10A
Rated Power	750W
Noise	≤65dB
External Size (W*D*H)	560*460*390mm
Net Weight (kg)	54
Package Size (W*D*H)	705*605*695mm
Gross Weight (kg)	70

ROTOR:

No.	Rotor	Capacity	Max. Speed (r/min)	Max. RCF (×g)
No. 1	Swing Rotor	4*250ml	4000	3500
	Adapter	4*12*5ml		