



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table Top Low Speed Refrigerated Centrifuge TLS-TL5R

FEATURES:

- Brushless DC motor
- LED/LCD display
- Stainless steel chamber
- Automatic electrical lid locking, over-speed, over-temperature and imbalance protection



MODEL SPECIFICATIONS:

Model	TLS-TL5R
Max. Speed	5000rpm
Max. RCF	4390xg
Max. Capacity	4*250ml
Speed Precision	±50rpm
Timing Range	1s~99h59min59s
Temp. Range	-20°C~40°C
Temp. Accuracy	±1°C
Refrigeration System	CFC-free refrigerant
Noise	≤65dB
Power Supply	AC110/220V, 50/60Hz
External Size (W*D*H)	640*680*380mm
Package Size (W*D*H)	780*800*520mm
Nett Weight (kg)	95
Gross Weight (kg)	120

ROTOR:

No.	Rotor	Capacity	Max. Speed (rpm)	Max. RCF (xg)
No. 1	Angle Rotor	12*10ml	5000	4390
No. 2	Swing Bucket Rotor	4*100ml	5000	4390
		4*2*100ml	4000	3500
		4*4*10ml		
		4*6*10ml		
		4*8*10ml		
		4*4*15ml		
		4*6*15ml		
		4*8*15ml		
		4*50ml	5000	4390
		4*2*50ml	4000	3500
No. 3	Swing Rotor	4*250ml	4000	3500
	Adapter	4*8*5ml		
		4*8*7ml		
		4*8*10ml		
		4*2*50ml		
		4*100ml		
No. 4	Swing Rotor	4*12*7/5ml	4000	2800
No. 5		4*18*7/5ml	4000	2800
No. 6		4*24*7/5ml	4000	2800
No. 7		Microplate	2*2*48 Well	4000
No. 8	2*2*96 Well		4000	2300