



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Urine Analyzer MA-100U

FEATURES:

- LCD display, keyboard operation
- Save data automatically
- Date and number can be entered
- The colour and limpiness can be entered
- Single- and continuous tests are available
- Adopt luminescence source with high brightness to minimize interference of natural light



MODEL SPECIFICATIONS:

Model	MA-100U
Throughput	120 tests/hr
Display	LCD display
Test Items	pH, Nitrite, Protein, Glucose, Occult Blood, Urobilinogen (URO), Bilirubin (BIL), Ketone (KET), Leukocytes (LEU), Glucose (GLU), Specific Gravity (SG), Ascorbic Acid, Color, Limpidness
Principle	Reflection of high-brightness cold light
Printer	Built-in thermal printer, paper type 57mm
Storage	1 000 patient data
Working Environment	Temperature: 0~40°C, humidity: 10~80%
Interface	Standard RS-232 port
Power Supply	AC220V±10% 50/60Hz; 110V±10% 60Hz
Package Size (W*D*H)	470*410*230mm
Gross Weight (kg)	3