



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Urine Analyzer MA-200U / MA-300U

FEATURES:

- Convenient, fast, accurate, simple
- Built-in thermal printer
- RS-232 interface can be connected to other systems, eg urine sediment system
- Fault diagnosis function
- Traceability: each inspection result is photographed and archived
- Environmental temperature monitoring system reduces the influence of the laboratory environment on the test results (only MA-300U)





MA-200U



MA-300U

MODEL SPECIFICATIONS:

Model	MA-200U	MA-300U
Measurement Principle	Albedometry	Artificial Intelligence + Computer Perception Technology
Test Items	14 items, including GLU, BIL, KET, SG, pH, BLD, PRO, URO, NIT, LEU, VC, CRE, CAL, MAL	
Test Speed	60-120 tests/h	514 tests/h
Display	LCD	7-inch touch screen
Wavelength	525-660nm	/
Working Environment	Temperature: 15~35°C, humidity: <85%	
Interface	RS-232 interface	
Storage	2 000 patient data	Mass storage
Printer	Built-in thermal printer	
Leakage Proof	Available, avoid cross contamination	
Power Supply	AC110V/220V±10% 50/60Hz	100~240V±10%, 50/60Hz
Package Size (W*D*H)	500*380*250mm	430*420*380mm
Gross Weight (kg)	4.5	5.5