

WATER PURIFIERS



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Water Purifier Consumables:

Model Type	PP Filter	5µm PP Deep Filter	20" PP Filter	PP Cotton Filter	Active Carbon Filter	Granular Active Carbon Filter	Active Carbon Block Filter	12L Deionized Purification Column (imported resin)	Deionized Purification Column	12L Ultra-Purified Column (imported resin)	Ultra-Purified Column (imported resin)	100/200/300/400/600 GPD RO Membrane	Organic Matter Purification Column	0.2µm Imported PES Terminal Filter	MWCO 5000D UF Ultra-Filtration Module	Dual Wavelength (185nm & 254nm) UV Lamp
Replacement term/set (months)	2-6	5-6	2-3	5-6	4-6	2-6	2-6	About 20000L pure water/set	About 1000L pure water/set	About 20000L pure water/set	About 1000L pure water/set	12-24	About 9000L pure water/set	/	/	About 9000 hours
MA-II(15L/H)	1	/	/	/	2	/	/	/	/	/	/	100*1	/	/	/	/
MA-II(30L/H)	/	/	/	1	2	/	/	/	/	/	/	200*1	/	/	/	/
MA-III (15L/H)	1	/	/	/	2	/	/	/	2	/	/	100*1	/	/	/	/
MA-III (30L/H)	/	1	/	/	2	/	/	/	3	/	/	200*1	/	/	/	/
MA-IV (45L/H)	/	/	1	/	/	1	1	/	/	/	/	300*1	/	/	/	/
MA-IV (63L/H)	/	/	1	/	/	1	1	/	/	/	/	400*1	/	/	/	/
MA-IV (94L/H)	/	/	1	/	/	1	1	/	/	/	/	600*1	/	/	/	/
MA-IV (125L/H)	/	/	1	/	/	1	1	/	/	/	/	400*2	/	/	/	/
MA-V (45L/H)	/	/	1	/	/	1	1	1	/	/	/	300*1	1	1	/	/
MA-V (63L/H)	/	/	1	/	/	1	1	1	/	/	/	400*1	1	1	/	/
MA-V (94L/H)	/	/	1	/	/	1	1	1	/	/	/	600*1	1	1	/	/

MA-V (125L/H)	/	/	1	/	/	1	1	1	/	/	/	400*2	1	1	/	/
MA-VI (45L/H)	/	/	1	/	/	1	1	/	/	1	2	300*1	1	1	/	/
MA-VI (63L/H)	/	/	1	/	/	1	1	/	/	1	2	400*1	1	1	/	/
MA-VI (94L/H)	/	/	1	/	/	1	1	/	/	1	2	600*1	1	1	/	/
MA-VI (125L/H)	/	/	1	/	/	1	1	/	/	1	2	400*2	1	1	/	/
MA-VII (45L/H)	/	/	1	/	/	1	1	/	/	1	2	300*1	1	1	1	/
MA-VII (63L/H)	/	/	1	/	/	1	1	/	/	1	2	400*1	1	1	1	/
MA-VII (94L/H)	/	/	1	/	/	1	1	/	/	1	2	600*1	1	1	1	/
MA-VII (125L/H)	/	/	1	/	/	1	1	/	/	1	2	400*2	1	1	1	/
MA-VIII (45L/H)	/	/	1	/	/	1	1	/	/	1	2	300*1	1	1	/	1
MA-VIII (63L/H)	/	/	1	/	/	1	1	/	/	1	2	400*1	1	1	/	1
MA-VIII (94L/H)	/	/	1	/	/	1	1	/	/	1	2	600*1	1	1	/	1
MA-VIII (125L/H)	/	/	1	/	/	1	1	/	/	1	2	400*2	1	1	/	1
MA-IX (45L/H)	/	/	1	/	/	1	1	/	/	1	2	300*1	1	1	1	1
MA-IX (63L/H)	/	/	1	/	/	1	1	/	/	1	2	400*1	1	1	1	1
MA-IX (94L/H)	/	/	1	/	/	1	1	/	/	1	2	600*1	1	1	1	1
MA-IX (125L/H)	/	/	1	/	/	1	1	/	/	1	2	400*2	1	1	1	1

Model Type	PP Filter	T33 Filter Element	Active Carbon Filter	Granular Active Carbon Filter	Precision Filter	20" Compressed Activated Carbon Filter	Capsule Filter	Deionized Purification Column	50/300/400/4021 GPD RO Membrane
Replacement term/set (months)	4-6	4-6	4-6	4-6	4-6	4-6	4-6	/	12
MA-I-10L	10"*1	/	2	/	/	/	/	3000L*2	50*1
MA-I-15L	10"*1	/	2	/	/	/	/	3000L*2	50*1
MA-II-20L	10"*1	/	1	/	/	/	/	5000L*1	300*1
MA-II-30L	10"*1	/	1	/	/	/	/	5000L*1	300*1
MA-II-40L	10"*1	/	1	/	/	/	/	5000L*1	300*1
MA-II-60L	/	/	/	/	1	/	/	8000L*1	300*2
MA-II-80L	/	/	/	/	1	/	/	8000L*1	300*2
MA-II-100L	/	/	/	/	1	/	/	8000L*1	400*2
MA-II-120L	20"*1	/	/	1	/	1	/	20000L*2	400*2
MA-II-150L	20"*1	/	/	1	/	1	/	20000L*2	400*2
MA-II-200L	20"*1	/	/	1	/	1	/	20000L*2	400*2
MA-V-45L	20"*1	1	/	1	/	1	1	20000L*1	300*1