

OVENS



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Constant Temperature Drying Oven MA-T30C / MA-T70C / MA-T140C / MA-T270C

FEATURES:

- External material: cold rolled steel with anti-bacteria powder coating
- Stainless steel inner chamber, with adjustable shelves
- PID control with LED display, soft touch operational key, timing function and over-temperature protection
- Door with double layer tempered glass observation window, adjustable airtight buckle lock and silicone rubber seal



MODEL SPECIFICATIONS:

Model	MA-T30C	MA-T70C	MA-T140C	MA-T270C
Capacity	43L	70L	137L	270L
Temp. Range	50~200°C			
Temp. Fluctuation	±1°C			
Ambient Temp.	5~40°C			
Timing Range	1~9999min			
Shelves No.	2 pcs		3 pcs	
Consumption	850W	1000W	2000W	3000W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)			
Internal Size (W*D*H)	350*350*350mm	400*390*450mm	500*500*550mm	600*600*750mm
External Size (W*D*H)	635*615*552mm	680*653*631mm	830*665*772mm	930*765*972mm
Package Size (W*D*H)	780*750*700mm	870*700*800mm	960*780*950mm	1060*880*1150mm
Gross Weight (kg)	60	70	105	145