



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

**Constant Temperature Water Tank**

**FEATURES:**

- Microprocessor control, LED display
- Timing function and over-temperature protection
- External material: cold rolled steel with anti-bacteria powder coating
- Stainless steel inner chamber



**MODEL SPECIFICATIONS:**

Model	MA-WT42	MA-WT60
Temp. Range	RT+5~100°C	
Temp. Fluctuation	≤±0.5°C	
Consumption	800W	1500W
Power Supply	AC110/200V±10%, 50/60Hz	
Internal Size (L*W*H) mm	420*180*180	600*300*180
External Size (L*W*H) mm	580*250*270	760*370*270
Package Size (L*W*H) mm	640*310*310	820*440*310
Net Weight (kg)	10	15
Gross Weight (kg)	11	17