



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

**Drying Oven / Incubator (Dual Use) MA-D30 / MA-D50 / MA-D70 / MA-D140 / MA-D240**

**FEATURES:**

- Over-temperature protection
- Microprocessor temperature control with LED display
- Stainless steel inner chamber, adjustable shelves
- Leakage protector





MA-D50



MA-D30



MA-D70 / MA-D140 / MA-D240

**MODEL SPECIFICATIONS:**

Model	MA-D30	MA-D50	MA-D70	MA-D140	MA-D240
Capacity	30L	50L	70L	140L	240L
Temp. Range	Incubator: RT +5~80°C; Drying Oven 80~200°C				
Temp. Precision	0.1°C				
Temp. Fluctuation	±0.5°C				
Ambient Temp.	5~40°C				
Timing Range	1~9999min				
Shelves No.	2 pcs				
Consumption	550W	850W	1100W	1500W	2050W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)				
Internal Size (W*D*H)	340*325*290mm	420*350*350mm	400*375*470mm	450*557*520mm	500*615*720mm
External Size (W*D*H)	490*520*690mm	710*520*520mm	540*550*870mm	600*740*920mm	650*800*1130mm
Package Size (W*D*H)	620*630*830mm	810*640*600mm	670*660*1040mm	730*860*1080mm	780*920*1270mm
Gross Weight (kg)	55	50	75	85	110