



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Forced Air Drying Oven MA-T25F / MA-T50F / MA-T105F / MA-T200F

#### FEATURES:

- External material: cold rolled steel with anti-bacteria powder coating
- Stainless steel inner chamber, with adjustable shelves
- PID control with LED display
- Double layered glass observation window, adjustable airtight buckle lock
- Over-temperature protection



**MODEL SPECIFICATIONS:**

Model	MA-T25F	MA-T50F	MA-T105F	MA-T200F
Capacity	25L	54L	108L	212L
Temp. Range	50~200°C			
Temp. Precision	0.1°C			
Temp. Fluctuation	±1°C			
Ambient Temp.	5~40°C			
Timing Range	1~9999min			
Shelves No.	2 pcs			
Consumption	650W	1000W	1600W	2500W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)			
Internal Size (W*D*H)	298*290*274mm	419*332*348mm	550*348*550mm	602*498*602mm
External Size (W*D*H)	580*554*455mm	708*595*532mm	830*620*735mm	880*775*785mm
Package Size (W*D*H)	710*660*620mm	840*690*700mm	960*720*900mm	1010*880*940mm
Gross Weight (kg)	51	66	115	145