



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

High Temperature Drying Oven MA-H50F / MA-H90F / MA-H216F

FEATURES:

- Mirror stainless steel inner chamber
- Optional: Test hole diameter 50mm on the left side
- Optional: Independent temperature protector
- Optional: RS485 interface for software connection



MODEL SPECIFICATIONS:

Model	MA-H50F	MA-H90F	MA-H216F
Capacity	50L	100L	2006L
Temp. Range	50~400°C		
Temp. Precision	0.1°C		
Temp. Fluctuation	±0.5°C		
Ambient Temp.	5~40°C		
Timing Range	1~9999min		
Shelves No.	2 pcs		
Consumption	2500W	3000W	3500W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)		380V, 50Hz
Internal Size (W*D*H)	355*380*400mm	450*450*480mm	600*600*610mm
External Size (W*D*H)	810*640*640mm	900*710*710mm	1060*860*860mm
Package Size (W*D*H)	940*800*800mm	1030*870*870mm	1190*1020*1020mm
Gross Weight (kg)	70	110	150