



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Oil Bath

#### FEATURES:

- LED display
- Automatic protection for over-heating
- External: steel plate with baking varnish; Internal: imported stainless steel
- Advanced recycle pump system ensures evenness of temperature (MA-OB2)



#### MODEL SPECIFICATIONS:

Model	MA-OB1	MA-OB2
Temp. Range	RT~300°C	
Temp. Fluctuation	≤1°C	
Temp. Control Precision	≤1°C	
Heating Power	2000W	
Power Supply	110/220V±10%, 50/60Hz	
Circulatory System	/	Yes
Circulation Pump Flow Rate	/	6L/min
Workroom Cubage	280*280*300mm <sup>3</sup>	300*300*280mm <sup>3</sup>
Package Size (W*D*H) mm	520*420*400	710*460*570
Gross Weight (kg)	15	33