



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Thermostatic Shaking Water Bath

FEATURES:

- Microprocessor control, LED display
- High precision of controlling temperature
- Advanced multipurpose stainless steel spring clip



MA-SWB A/C



MA-SWB B

MODEL SPECIFICATIONS:

Model	MA-SWB A	MA-SWB B	MA-SWB C
Temp. Range	RT~100°C		RT~99.9°C
Oscillation Amplitude	30mm	24mm	20mm
Rotary Speed	0~320rpm	50~280rpm	0~300rpm
Shaking Mode	Rotation oscillation, continuous	Rotation oscillation & reciprocating, continuous	Reciprocating, continuous
Timing	0~120min		
Consumption	1500W	1550W	1500W
Power Supply	AC110/200V±10%, 50/60Hz		
Internal Size (L*W*H) mm	490*390*170	490*390*150	490*390*170
External Size (L*W*H) mm	720*480*520		
Package Size (L*W*H) mm	820*596*560	800*590*580	800*590*560
Gross Weight (kg)	48	52	48