



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Thermostatic Water Bath

FEATURES:

- Stainless steel inner chamber and upper cover
- PID microprocessor controller, LED display
- With universal fixture for MA-1L4H and MA-1L6H



MA-1L1H



MA-2L4H



MA-2L8H

MODEL SPECIFICATIONS:

Model	MA-1L1H	MA-1L2H	MA-1L4H	MA-1L6H	MA-2L4H	MA-2L6H	MA-2L8H
Capacity	3L	6L	13L	20L	11L	18L	23L
Temp. Range	RT~100°C						
Temp. Fluctuation	±0.5°C						
Temp. Accuracy	≤±1°C						
Time Range	0~999min						
Power	300W	500W	1000W	1500W	1000W	1500W	2000W
Power Supply	AC110/200V±10%, 50/60Hz						
Internal Size (mm)	160*160*130	300*160*130	600*180*130	900*180*130	300*300*130	470*300*130	600*300*130
External Size (mm)	200*200*220	345*210*220	770*230*170	1060*230*170	345*345*220	500*340*220	770*370*170
Package Size (mm)	260*260*280	460*260*280	830*300*250	1150*300*250	400*400*280	600*400*280	830*460*250
Net Weight (kg)	2	4	10	13	6.5	9	12
Gross Weight (kg)	3	5	11	14.5	7.5	10	13