



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Constant-Temperature Incubator MA-H30 / MA-H50 / MA-H80 / MA-H160 / MA-H270

#### FEATURES:

- Intelligent temperature controller, based on a micro-processor with PID parameters and LCD display
- Double door design to minimise temperature fluctuations
- Toughened glass inner door to observe chamber conditions without affecting chamber temperature
- Optional accessory: independent temperature protector



MA-H30



MA-H50/MA-H80/MA-H160/MA-H270

**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>MA-H30</b>	<b>MA-H50</b>	<b>MA-H80</b>	<b>MA-H160</b>	<b>MA-H270</b>
Capacity	30L	50L	80L	160L	270L
Temp. Range	RT+5~65°C				
Temp. Fluctuation	±0.3°C				
Time Setting	1~9999mins				
Alarm	Over-temperature alarm; audible and visual alarm				
Shelves	2 pcs/adjustable			3 pcs/adjustable	
Consumption	200W	210W	280W	380W	570W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110/220V±10%, 60Hz (optional)				
Internal Size (W*D*H) mm	320*300*315	350*350*410	400*400*500	500*500*650	600*600*750
External Size (W*D*H) mm	650*470*510	510*520*800	580*570*890	680*670*1050	770*780*1150
Package Size (W*D*H) mm	780*600*680	600*610*850	670*650*940	800*790*1200	890*890*1300
Gross Weight (kg)	45	55	57	90	130