

INCUBATORS



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za


Constant Temperature Incubator MA-H30II / MA-H48II / MA-H64II / MA-H123II / MA-H230II

FEATURES:

- Large LED display, shows multiple groups of data at a time
- Timing and over-temperature protection
- Double layered glass observation window



MA-H30II/MA-H48II/MA-H64II/MA-H123II



MA-H230II

MODEL SPECIFICATIONS:

Model	MA-H30II	MA-H48II	MA-H64II	MA-H123II	MA-H230II
Capacity	30L	45L	65L	125L	230L
Temp. Range	RT+5~80°C				
Temp. Fluctuation	±0.2°C				
Time Setting	1~9999mins				
Shelves	2 pcs / adjustable				
Consumption	150W	250W	500W	600W	700W
Power Supply	AC220V±10%, 50/60Hz (standard), AC110V±10%, 60Hz (optional)				
Internal Size (W*D*H) mm	310*310*310	350*350*350	400*360*450	500*450*550	600*500*750
External Size (W*D*H) mm	460*510*695	500*550*735	550*550*840	636*680*915	730*670*1220
Package Size (W*D*H) mm	550*570*765	620*585*800	640*640*905	730*720*1000	900*800*1350
Gross Weight (kg)	51	56	62	71	105