

## INCUBATORS



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Constant Temperature Incubator MA-H50II / MA-H80II / MA-H160II / MA-H270II

#### FEATURES:

- LCD display
- Micro-computer temperature controller
- Double door design and magnetic door seal
- Incubator equipped with castors for easy movement

#### SAFETY FUNCTIONS:

- Over-temperature protection, alarm for temperature sensor malfunction
- Parameters can be saved automatically during power outages or system errors



MA-H50II-03

**MODEL SPECIFICATIONS:**

Model	MA-H50II	MA-H80II	MA-H160II	MA-H270II
Capacity	48L	84L	160L	270L
Temp. Range	RT+5~60°C			
Temp. Fluctuation	±0.3°C			
Time Range	1~9999mins			
Shelves	2 pcs/adjustable		3 pcs/adjustable	
Consumption	437W	502W	580W	688W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)			
Internal Size (W*D*H) mm	350*348*400	410*405*510	505*455*700	600*600*750
External Size (W*D*H) mm	510*548*740	570*590*850	673*645*1033	768*795*1083
Package Size (W*D*H) mm	590*650*930	730*700*1030	800*800*1210	930*940*1250
Gross Weight (kg)	50	75	99.5	150