



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Lighting Incubator MA-L200BK

FEATURES:

- Large LCD screen
- Forced-air convection
- Micro-processor control
- Universal castors with lock for easy movement
- Double door design, outer door with viewing window
- Polished stainless steel chamber

SAFETY FUNCTIONS:

- Over-temperature alarm
- Independent limit temperature alarm system
- Leakage and over-voltage protector

OPTIONAL ACCESSORIES:

- RS485 interface
- Temperature controlled switch



Large LCD screen



Mechanical Lock



Test port



Outer door with viewing window

MODEL SPECIFICATIONS:

Model	MA-L200BK
Capacity	200L
Lighting Type	LED cold light source
Temp. Range	With illumination: 10~60°C; Without illumination: 5~60°C
Temp. Resolution	0.1°C
Temp. Fluctuation	±0.5°C
Temp. Deviation	±1.0°C
Lighting Illumination	0~12000 lux (5-grade adjustable)
Ambient Temp.	0~40°C
Continuous Working Time	No less than 180h
Refrigerant	R134a (CFC free)
Shelves	Standard: 2 pcs
Heating Power	900W
Cooling Power	300W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)
Internal Size (W*D*H) mm	548*472*778
External Size (W*D*H) mm	760*750*1500
Package Size (W*D*H) mm	883*870*1645
Gross Weight (kg)	200