

INCUBATORS



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Medicine Stability Test Chamber (Economic Type) MA-MS250 / MA-MS300 / MA-MS400

FEATURES:

- PID accurate, reliable temperature and humidity control
- Built-in fan forms air circulation, improving temperature uniformity inside the cabinet
- Built-in water tank and automatic water replenishment interface, automatically disconnects heater in the event of lack of humidification water
- Accessories:
 - Standard: Micro thermal printer
 - Optional: USB data storage, micro-needle printer, external temperature protection device



MA-MS250-03

MODEL SPECIFICATIONS:

Model	MA-MS250	MA-MS300	MA-MS400
Display	5.7-inch LCD		
Capacity	248L	290L	402L
Temp. Range	0~65°C		
Temp. Fluctuation	±0.5°C		
Temp. Uniformity	±2°C		
Humidity Control Range	40~95% RH (Depends on temperature beyond 10°C)		
Humidity Fluctuation	±3% RH		
Illumination	6000 Lux, 5 level adjustable, 1-side illumination		
Shelves	2 pcs/adjustable		
Timing Range	1~99min		
Refrigeration System	R134a refrigerant, CFC-free		
Alarm	High / Low temperature alarm, temperature and humidity probe damage alarm		
Consumption	2400W	2900W	3400W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)		
Chamber Size (W*D*H) mm	470*449*1178	490*464*1278	560*600*1198
External Size (W*D*H) mm	672*892*1825	691*907*1925	764*1032*1860
Package Size (W*D*H) mm	800*1010*1990	790*1040*2090	870*1140*2040
Gross Weight (kg)	158	178	200