



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Medicine Stability Test Chamber MA-MS120A / MA-MS288A / MA-MS500A

Provides a stable environment of temperature, humidity and illumination for the pharmaceutical industry medicine invalidation evaluation.

FEATURES:

- PID accurate, reliable temperature and humidity control
- Built-in fan forms air circulation, improving temperature uniformity inside the cabinet
- Built-in water tank and automatic water replenishment interface, automatically disconnects heater in the event of lack of humidification water
- 50 programmes for the pharmaceutical application
- 6/day historical data USB export function
- Optional: Printer, external temperature protection device



MA-MS120A-03

MODEL SPECIFICATIONS:

Model	MA-MS120A	MA-MS288A	MA-MS500A
Display	7-inch LCD		
Capacity	173L	288L	516L
Temp. Range	-20~65°C		
Temp. Fluctuation	±0.5°C		
Temp. Uniformity	±2°C		
Humidity Control Range	40~95% RH		
Humidity Fluctuation	±3% RH		
Illumination	4500 ± 500 Lux, 4 level adjustable, 1-side illumination		
Shelves	2 pcs/adjustable		
Timing Range	1~9999min		
Refrigeration System	R404a refrigerant, CFC-free		
Alarm	High / Low temperature alarm, temperature and humidity probe damage alarm		
Consumption	2100W	2600W	3500W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)		
Chamber Size (W*D*H) mm	550*450*700	600*600*800	800*710*910
External Size (W*D*H) mm	700*842*1585	800*992*1685	1000*1072*1795
Package Size (W*D*H) mm	810*1080*1780	950*1160*1940	1170*1460*2050
Gross Weight (kg)	138	168	250