



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Touch Screen Constant Temperature Incubator MA-H54BK (D/G) / MA-H88BK (D/G) / MA-H160BK (D/G)

FEATURES:

- Large LCD touch screen
- D-type: with inner glass door; G-type: with viewing window
- Micro-processor controller
- Stainless steel chamber

SAFETY FUNCTIONS:

- Over-temperature alarm
- Independent temperature limit alarm system
- Leaking and over-voltage protector



Double layered door
(D-type)



Glass Viewing Window Door
(G-type)

Standard Accessories:



Test Port

Optional Accessories:



Base Stand



Universal Casters



UV Lamp

MODEL SPECIFICATIONS:

Model	MA-H54BK (D/G)	MA-H88BK (D/G)	MA-H160BK (D/G)
Capacity	54L	88L	160L
Temp. Control Method	PT100		
Temp. Range	RT+5~85°C		
Temp. Stability	±0.5°C (at 37°C)		
Temp. Uniformity	±1°C (at 37°C)		
Ambient Temp.	5~30°C, recommended 25°C ±2°C		
Timing Range	0~999h		
Shelves	2 pcs/adjustable		
Consumption	345W		495W
Power Supply	AC220V±10%, 50/60Hz (standard); AV110V±10%, 60Hz (optional)		
Internal Size (W*D*H) mm	395*346*395	500*364*500	558*518*600
External Size (W*D*H) mm	D-type: 536*588*681 G-type: 522*626*587	D-type: 650*573*793 G-type: 627*621*692	D-type: 718*732*890 G-type: 685*779*790
Package Size (W*D*H) mm	D-type: 606*648*818 G-type: 592*676*724	D-type: 720*630*930 G-type: 697*671*829	D-type: 788*792*1027 G-type: 755*829*927
Gross Weight (kg)	D-type: 51 G-type: 44	D-type: 65 G-type: 51	D-type: 89 G-type: 69