

## HORIZONTAL AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Horizontal Autoclave BKQ-B100H / BKQ-B150H / BKQ-B200H / BKQ-B300H

#### APPLICATION RANGE:

For the sterilization of unpacked sold instruments in stomatology, ophthalmology, operation rooms and CSSD.

#### FEATURES:

- Sterilizing programs for fabric, instruments, rubber and liquid
- Touch screen displays temperature, pressure, time, running status and failure warning
- Post-drying function makes the residual humidity of the sterilized articles less than 3%
- Equipped with printer to print date, time and process parameters
- Built-in efficient steam generator, microcomputer temperature control
- Double door optional



BKQ-B100H-03

**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>BKQ-B100H</b>	<b>BKQ-B150H</b>	<b>BKQ-B200H</b>	<b>BKQ-B300H</b>
Capacity	100L	135L	185L	300L
Chamber Size (mm)	φ 400*850	φ 500*710	φ 500*950	φ 632*1000
Chamber Material	SUS304			
Max. Working Pressure	0.23Mpa			
Working Temp.	115~138°C			
Temp. Accuracy	±0.1°C			
Consumption	10KW	11KW		17KW
Power Supply	380V, 50Hz			
External Size (W*D*H) mm	1250*730*1660	1110*750*1742	1350*750*1742	1400*890*1780
Package Size (W*D*H) mm	1500*900*2050	1370*1000*2060	1510*890*2060	1700*1050*2090
Gross Weight (kg)	310	470	540	830