

HORIZONTAL AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Horizontal Cylindrical Pressure Autoclave BKQ-H150 / BKQ-H200 / BKQ-H300 / BKQ-H400 / BKQ-H500

FEATURES:

- Drying function
- Interlock function
- Automatic sterilization process
- 304 stainless steel material chamber
- Safety valve for over pressure protection
- Over-temperature, over-pressure protection device
- Water level sensor, power will be cut-off
- Double door optional



BKM-H150-03

MODEL SPECIFICATIONS:

Model	BKQ-H150	BKQ-H200	BKQ-H300	BKQ-H400	BKQ-H500
Capacity	150L	200L	280L	400L	500L
Chamber Size (mm)	φ 440*1000	φ 515*1000	φ 600*1000	φ 700*1100	φ 700*1300
Working Pressure	0.22MPa				
Working Temp.	134°C				
Temp. Range	105~134°C				
Sterilization Time	0~60min				
Drying Time	0~60min				
Temp. Display Accuracy	≤±2°C				
Consumption	9KW		12KW	18KW	
Power Supply	220V/380V, 50/60Hz		380V, 50Hz		
External Size (W*D*H) mm	1400*600*1300	1400*670*1650	1400*770*1780	1430*880*1830	1800*900*1820
Package Size (W*D*H) mm	1550*750*1850	1560*820*1850	1680*920*2100	1600*1050*2100	1850*1050*2100
Gross Weight (kg)	320	350	465	530	580