

HORIZONTAL AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Horizontal Pulse Vacuum Autoclave BKQ-Z100(H) / BKQ-Z150(H) / BKQ-Z200(H) / BKQ-Z300(H)

APPLICATION RANGE:

Can be widely used in stomatology, ophthalmology, operation rooms and other medical institutions. It is suitable for all wrapped or unwrapped solid instruments, A-class cavity instruments (dental hand pieces and endoscopes), implantable instruments, dressing fabric and rubber tubes etc.

FEATURES:

- Built-in steam generator, microcomputer temperature control
- Equipped with printer to print date, time and process parameters
- Four-level authority management system
- Sterilization programs for fabric, instruments, rubber, liquid, cavity, B&D tests, PCD tests and leakage tests
- Intelligent drying system to guarantee the residual humidity of instrument and dressing no more than 0.2% and 1% respectively
- Overheat auto protection device, dual door safety lock, over pressure relief valve, short circuit protection
- LCD screen displays temperature, pressure, time, running status, failure and warning codes
- Double door optional



BKQ-Z100(H)-03

MODEL SPECIFICATIONS:

Model	BKQ-Z100(H)	BKQ-Z150(H)	BKQ-Z200(H)	BKQ-Z300(H)
Capacity	100L	135L	185L	300L
Chamber Size (mm)	φ 400*850	φ 500*710	φ 500*950	φ 632*1000
Chamber Material	SUS304			
Max. Working Pressure	0.23Mpa			
Vacuum Limit	-0.08Mpa			
Working Temp.	115~138°C			
Temp. Accuracy	±0.1°C			
Consumption	11KW	11.5KW		18KW
Power Supply	380V, 50Hz			
External Size (W*D*H) mm	1250*730*1660	1110*750*1742	1350*750*1742	1400*890*1780
Package Size (W*D*H) mm	1500*900*2050	1370*1000*2060	1510*890*2060	1700*1050*2090
Gross Weight (kg)	310	470	540	830