

TABLE TOP AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table-Top Autoclave Class B-Series



BKM-Z18BIII/24BIII/45BIII



BKM-Z60BIII/80BIII

FEATURES:

- High pressure sterilizer, using steam as its sterilization medium, which is fast, safe and economic. Commonly used in stomatology, ophthalmology, operating rooms and CSSD for sterilization of wrapped or unwrapped instruments, fabric, utensils, culture medium, unsealed liquid etc
- Built-in open type water tank – water quality monitor, water filter inside, water level sensor, a fully injected water tank supports repeated program running
- LCD screen displays temperature, pressure, time, operating status, error and alarm information

- Multiple security protective device – over-temperature automatic protective device, multiple control and protection for steam generator, door safety interlock protection and dual chamber over-pressure protection, electronic circuit safety device
- Automatic door – BKM-Z18BIII/24BIII/45BIII series has automatic door structure, the BKM-60BIII/80BIII series has multipoint synchronous linkage pressing technology
- Efficient vacuum, low working noise, vacuum lowest limit can reach over -90KPa, vacuum in high speed, thorough sterilization
- Multiple program types – sterilization programs: unwrapped, wrapped, rubber, user 01, N-quick, B-quick, prion test program, B&D Helix, leakage test, auxiliary program: preheat, drying, washing
- High efficiency steam generator – BKM-Z18BIII/24BIII/45BIII series has built-in instant steam generator, BKM-Z60BIII/80BIII series has built-in energy storing type steam generator, producing steam in high speed and large quantity
- Other features: Pulse vacuum; Heating; Sterilization; Exhausting; Vacuum drying totally automatic; Chamber is heated by film heater; Record type: built-in printer (standard); Electronic storage: USB flash disk (optional)

Loading tray:



Type A



Type B



Type C (optional for 45L/60L/80L)

MODEL SPECIFICATIONS:

| Model | BKM-Z18BIII | BKM-Z24BIII | BKM-Z45BIII | BKM-Z60BIII | BKM-Z80BIII |
|--------------------------|-------------------------|-------------|---------------------|--------------|--------------|
| Capacity | 18L | 24L | 45L | 60L | 80L |
| Chamber Size (mm) | φ 250*350 | φ 250*450 | φ 316*618 | φ 386*595 | φ 386*717 |
| Design Pressure | -0.1/0.3MPa | | | | |
| Design Temp. | 144°C | | | | |
| Temp. Range | 105~138°C | | | | |
| Chamber Material | SUS304 | | | | |
| Standard Loading Tray | Type A | | Type B | | |
| Consumption | 2.6KW | | 5.8KW | 7.4KW | |
| Power Supply | AC110/220V±10%, 50/60Hz | | AC220V±10%, 50/60Hz | | |
| External Size (W*D*H) mm | 530*560*450 | 530*660*450 | 640*830*550 | 730*1002*560 | 730*1002*560 |
| Package Size (W*D*H) mm | 720*780*680 | 720*880*700 | 750*1015*730 | 830*1180*795 | 830*1180*795 |
| Gross Weight (kg) | 70 | 80 | 140 | 185 | 205 |