

TABLE TOP AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table-Top Autoclave Class B-Series

FEATURES:

- High pressure sterilizer, using steam as its sterilization medium, which is fast, safe and economic. Commonly used in medical health care institutions or R&D institutions for sterilization of wrapped or unwrapped instruments, fabric, utensils, culture medium, unsealed liquid etc
- LCD touch screen displays temperature, pressure, time, operating status, error- and alarm information. There are different colour changes in each stage of sterilization
- Multiple safety protection devices – automatic over-temperature protection device, multiple control and protection of evaporator, door safety interlock protection, double over-pressure protection etc
- Mobile temperature sensor
- Automatic door structure – sterilizer adopts automatic door structure, fast and convenient to operate
- Vacuuming system – high efficiency vacuum extraction system, low working noise, over -90kPa extreme vacuum, fast extraction speed, more complete sterilization
- Variety of programs –exposed procedures, packaging procedures, rubber procedures, dressing procedures, liquid procedures, liquid customization, solid customization, medium sterilization procedures, B&D testing, vacuum testing and media dissolution procedures



BKM-Z45B-03

- Other features: Heating, sterilizing, automatic exhausting and drying; Chamber is heated by film heater; 0.22 μ m high efficiency air filter to avoid re-contamination of sterilized materials with filtration efficiency over 99.5%; Water filter in tank and chamber; Built-in stainless steel water tank, inner cycle of steam and water; Standard built-in printer, Standard USB port

MODEL SPECIFICATIONS:

Model	BKM-Z45B	BKM-Z60B	BKM-Z80B
Capacity	45L	60L	80L
Chamber Size (mm)	ϕ 316*621	ϕ 386*520	ϕ 386*692
Design Pressure	-0.1/0.3MPa		
Sterilization Temp.	105°C ~ 134°C		
Chamber Material	SUS304		
Ambient Temp.	5°C~40°C		
Consumption	5.8KW	6.2KW	
Power Supply	AC110/220V \pm 10%, 50/60Hz		
External Size (W*D*H) mm	610*1020*610	805*940*700	805*1140*700
Package Size (W*D*H) mm	750*1200*880	960*1130*940	960*1330*940
Gross Weight (kg)	165	215	230

Note: It is applicable to altitude below 3,500m. Parts need to be changed when altitude exceeds 2,000m