

TABLE TOP AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Table-Top Autoclave Class N-Series

FEATURES:

- For use in microbiology, medicine, podiatry, tattooing, body piercing, veterinary science, mycology, dentistry and prosthetics. Vary in size and function depends on the media to be sterilized
- Adopting European N standard, with thermal vacuum and drying functions, it can be used for A-type hollow - and solid instruments. Thermodynamic vacuum guarantees elimination of cold air from the cavity, ensuring optimal effect of sterilization
- Specifically designed for the dental clinic field, it could sterilize 12 pieces simultaneously
- Microprocessor control, LED display and modularized panel
- Automatic malfunction detection system



BKM-Z12N



BKM-Z18N/24N

BKM-Z12N-04

MODEL SPECIFICATIONS:

Model	BKM-Z12N	BKM-Z18N	BKM-Z24N
Capacity	12L	18L	24L
Chamber Size (mm)	φ 200*360	φ 245*350	φ 245*450
Sterilization Class	Class N (according to EN 13060 standard)		
Sterilization Temp.	121°C, 134°C		
Drying Procedure	Heating dry		
Display	LED display		
Control Precision	Temperature: 1°C		
	Pressure: 0.1 bar		
Safety System	Hand lock door		
	Relief valve in case of over-pressure		
	Pressure- or temperature overload protection		
	Alarm for system failure, stop reminder		
Tray Holder	3 pcs SS trays on SS shelf		
Chamber	SUS304		
	Max. working pressure: 2.3 bar		
	Min. working pressure: 1.1 bar		
	Max. temperature: 145°C		
Ambient Temp.	5~40°C		
Noise	<50dB		
Consumption	1400W	1500W	1700W
Power Supply	AC110/220V±10%, 50/60Hz		
External Size (W*D*H) mm	560*450*400	560*450*400	670*450*400
Package Size (W*D*H) mm	660*525*495	660*525*495	770*525*495
Gross Weight (kg)	50	53	60