

## TABLE TOP AUTOCLAVE



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Table-Top Dental Sterilizer

#### FEATURES:

- Microcomputer control, LCD display and simple operation interface
- New designed type with opening water tank on top
- Small capacity design more suitable for clinic requirements
- Optional Accessory: mini printer



BKM-Z8B-03

**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>BKM-Z8B</b>	<b>BKM-Z12B</b>
Capacity	8L	12L
Chamber Size (mm)	φ 170*320	φ 192*320
Sterilization Class	Class B	
Sterilization Temp.	121°C, 134°C	
Chamber Material	SUS304	
Max. Working Pressure	0.25Mpa	
Max. Working Temp.	145°C	
Consumption	1.4KW	
Power Supply	AC110/220V±10%, 50/60Hz	
External Size (W*D*H) mm	540*370*360	
Package Size (W*D*H) mm	615*440*425	
Gross Weight (kg)	45	50