

## VERTICAL AUTOCLAVE



+27 60 999 9140


info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### Vertical Autoclave



BKQ-B50L-03

**FEATURES:**

- Especially designed for liquid and food sterilization in laboratories, pharmaceutical and food industry
- Standard built-in printer
- Deionized water interface to supply the steam generator
- F0 valued displayed and printed
- Rapid cooling water system for efficient sterilization of liquids
- Quick maintenance window
- Microcomputer control, fully automatic control water injection, temperature rise, sterilization etc
- LCD screen shows pressure, temperature, time, running status, error codes and operation curve etc
- Mobile probe detects the internal temperature of liquid directly, ensuring the sterilization effect of liquid program

**MODEL SPECIFICATIONS:**

Model	BKQ-B50L	BKQ-B75L
Capacity	50L	75L
Chamber Size (mm)	φ 386*490	φ 386*670
Chamber Material	SUS304	
Designed Pressure	0.28Mpa	
Designed Temp.	150°C	
Working Pressure	0.22Mpa	
Working Temp.	80°C~136°C	
Temp. Display Accuracy	0.1°C	
Noise	≤65dB	
Consumption	5.5KW	
Power Supply	AC110/220V±10%, 50/60Hz	
SS Basket	1 pc	
External Size (W*D*H) mm	700*610*1100	
Package Size (W*D*H) mm	800*715*1270	
Gross Weight (kg)	140	147