

CENTRIFUGES

MINI-CENTRIFUGES



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Mini-Centrifuge MA Mini-10K / MA Mini-10k+

Used for serum extraction from whole blood, extraction from different samples, micro blood cell separation, bio-sample treatment, samples quick sedimentation and PCR sample treatment.

FEATURES:

- Timing function
- Rotor with 8 positions
- Transparent cover
- Auto-stop when cover opens during running and will restart running after cover is closed
- Low noise, stable running
- Unique rotor clip design



MODEL SPECIFICATIONS:

Model	MA Mini-10K	MA Mini-10K+
Max. Speed	10000rpm	3000-10000rpm
Max. RCF	5400xg	500-5400xg
Rotor	8*2.0ml/1.5ml/0.5ml/0.2ml (0.2ml and 0.5ml adaptor) 2*8*0.2ml PCR centrifugal tube	
Noise	≤55dB	
Power Supply	AC110/220V±10%, 50/60Hz	
External Size (W*D*H)	180*160*121mm	
Package Size (W*D*H)	280*220*200mm	
Net / Gross Weight (kg)	1.4 / 3	